

START

START SESSION

Tap **START SESSION**.

START CALIBRATION

Tap **START CALIBRATION** and stand still with arms folded.

START SESSION

When calibration is complete, tap **START SESSION**.

ADJUST

Select Screen
to be Displayed



Workout session
information

Adjust Incline

Adjust Body weight

Adjust Speed

PAUSE | STOP



Tap **PAUSE** to pause the workout session. The screen will become grayed out and the **GO** button will appear.



Tap **STOP** to end the workout session.



Tap **GO** to resume the session.

EMERGENCY EXTRACTION PROTOCOLS

Before using the Anti-Gravity Treadmill, read the entire User Manual, which can be obtained at alterg.com or by contacting our Service Department at +1.510.270.5900.

In the unlikely event that a patient becomes immobilized in the Anti-Gravity Treadmill, please refer to the urgent steps below to secure the situation before calling +1510.270.5900 or referring to alterg.com or the User Manual for more detailed instructions and next steps.

IMPORTANT NOTES:

1. Do not unlock the cockpit until you read further.
2. The treadmill belt on the Via series does not lock in place, so secure it by holding your foot against it at the rear of the treadmill.

EMERGENCY SCENARIOS

There are two likely scenarios in which a patient may require urgent help with getting out of the Anti-Gravity Treadmill:

SCENARIO 1

The patient falls or becomes fatigued to the point of inability to exit the treadmill before the air chamber is inflated with air.

- Option 1:** If the patient can stand, stop the treadmill and secure the belt with your foot. Help the patient exit the treadmill.
- Option 2:** Prevent the belt from moving by securing it with your foot. Remove the bag rods at the treadmill base.
- Option 3:** Prevent the belt from moving by securing it with your foot. Cut the air chamber and help the patient exit the treadmill.

SCENARIO 2

The patient becomes fatigued to the point of inability to exit the treadmill after the air chamber is inflated.

Reduce the body weight percentage to less than 40% and have the patient sit down in place until they re-gain their strength and are able to exit on their own.